

## Dräger Perseus® A500 Breathing System Disassembly/Assembly



### Disassembly

1 Unscrew the CLIC adapter.



Screw the CLIC adapter.

2 Insert the cleaning disk into the adapter. Close the adapter so that it engages.



Remove the cleaning disk from the adapter. Close the adapter so that it engages.

3 Loosen and pull out the expiration port.



Push in the expiratory port and tighten it with the knurled nut.

4 Loosen and pull out the inspiratory port.



Push in the inspiratory port and tighten it with the knurled nut.

5 Remove the flow sensors.



Insert the flow sensors.

6 Turn the three levers counterclockwise about 120°.



Turn the three levers clockwise about 120°.

7 Lift off the upper part of the breathing system.



Place the upper part of the breathing system on the lower part.

8 Screw off the APL valve.



Align the APL valve correctly and tighten it with the knurled nut.

9 Remove the valves (red and blue).



Align the valves correctly. Insert the valves.

### Assembly

### Reprocessing instruction\*

#### Disinfection by immersion



Please use the disinfectants listed in the Instructions for use of the Dräger flow sensors.

#### Machine cleaning and disinfection



Position components so that all interior spaces and surfaces are completely flushed and water can drain off freely.

\* For machine cleaning and disinfection of the blower the adapter MK10136 is required.

#### Sterilization 134 °C (273.2 °F)



SpiroLife® flow sensors are not suitable for plasma or radiation sterilization.

#### Visual check



The following parts must be free of deposits:  
– Incident flow grills  
– Non-return valves (red and blue)  
If required, remove deposits on the valve plates of the non-return valves with a soft cloth. If the incident flow grills are damaged, contact DrägerService®.

\*This document is not a replacement for the Instruction for Use.