



Neonatal Ventilator | SLE2000



When the smallest thing matters

SLE2000: The Reliable Solution for Infant Ventilation



For many years, the SLE2000 has set the standard in neonatal ventilation.

Using a unique valveless system, the principle of operation of the SLE2000 helps to eliminate inadvertent PEEP, and aid in the total clearance of expired gases.

- Ability to maintain pressure waveforms at all rates
- Sensitive airway pressure trigger
- Built in oxygen analyser with continuous digital readout
- Excellent pressure square wave means rapid inflation
- High pressure dump on alarm
- Optional auxiliary blended flow outlet
- Choice of square or slow rise time (switch selectable)
- Can be used with Nitric Oxide therapy
- Low maintenance costs

The smallest thing does matter...

The SLE2000 is specifically designed for neonatal and infant patients. It allows the user to deliver conventional ventilation i.e. CPAP, CMV, PTV or SIMV.

Exceptional pneumatic performance is further enhanced by a patented valveless system.

The valveless system gives the SLE2000 many advantages over conventional ventilators including a less complicated patient circuit with no diaphragms or other obstructions in the exhalation limb of the circuit, helping to reduce expiratory resistance.

All of which makes it easier to assemble and clean and ultimately reduces the potential for errors or wrong assembly.

Customer Support

Throughout the 50 plus year history of SLE, the company has always put customers and patients first.

We believe that developing life-support equipment should not be taken lightly, and that only the highest quality engineering is acceptable for our products.

Downtime is not desirable in the modern NICU. Logical product design means that routine servicing can be carried out quickly and effectively.

Our teams of qualified engineers support our distributors and end-users world-wide to ensure that your SLE product is always available for optimum performance.

Clinical Education

We know that a ventilator is a vital part of your equipment and that you need to feel comfortable using it. That's why we employ experienced Clinical Specialists in our Clinical Education Department.

These Specialists support our world-wide network of distributors and can help you become more familiar with your SLE products.

Ordering details

Z2002/FL SLE2000 with single-use patient circuits (International)

Z0002 RS-232 serial port

Z0003 0-15 l/min blended output

Z0004/01 Factory fitted modification to allow BPM range of 126 - 250 BPM

Z0005 External moisture trap (on air supply only)

N2391 Patient Circuit-Reusable (10 mm) (Other types available on request)

N2587/000/050 Bacteria Filter-Single use (Box of 50)

N2029 Bacteria Filter - Autoclavable

N2186 Silencer - Autoclavable

N2190 Exhalation block - Autoclavable

N2002/FL Service manual

N2000/FL User manual

N2001/FL Quick Setup Guide

Compressor

A free-standing air compressor is also available for use with the SLE2000. Please see your local distributor for further details.

BC2188/400/15

Single use breathing circuit for use with SLE2000 and SLE2000 HFO infant ventilators. Temperature port 400 mm from ET manifold (single use). Circuit comes complete with water trap, filter connection kit and adaptors. (Box of 15)

BC6216

Nitric Oxide delivery kit, set of connectors (Paediatric delivery).

BC2508

Nebuliser kit (Paediatric delivery).

BC2110/KIT

Nitric Oxide adaptor kit for BC2188/400 breathing circuits (SLE2000 and SLE2000 HFO infant ventilators). (Box of 5)

N4110

Nitric Oxide dual hose scavenging filter assembly for SLE2000 infant ventilators. (Box of 1)

Ventilator Controls

Modes: Off/Alarm Test, CPAP, CMV, PTV, SIMV

Breathing Rates: 1-125 (126-250 with Option Z0004)

Max Inspiratory Time: 0.1-3.0 seconds (0.01-0.3 seconds with Option Z0004)

I.E. Ratio: 9.9:1 to 1:99 calculated from BPM and Ti settings

Pressure Wave: Slow or fast rise, switch selectable

Manual Breath: Gives single breath in CPAP, CMV, PTV, SIMV modes

PTV Sensitivity: Adjustable threshold

Oxygen Blender: 21-100% O₂ ± 3%

Pressure Display Switch: Max - Mean - Min

CPAP Regulator: Range: 0-15 mbar

Inspiratory Pressure Regulator: 0-60 mbar

Alarm Mute: 60 seconds

Alarm Reset: Resets all alarms (except system fail)

Indicator LEDs

Power: Green LED, Power on

System Fail: Red LED, Main processor system fail

Trigger Back up: Green LED, Indicates a machine delivered breath

Proximal airway LED

Pressure meter: Range -4 to +62 mbar

Digital

BPM Display: Breaths per minute

Insp Time Display: Inspiratory time

IE Ratio: Inspiration/ Expiration ratio

FiO₂: Oxygen percentage 21 - 100%

Pressure: 0-62 mbar

Gauge

Driving Pressure Gauge: 0-60 mbar plus PEEP level

Alarms

CPAP: Visual and intermittent audible

PIP / Cycle Fail / Low: Visual and intermittent audible

Fresh Gas Block or Leak: Visual and continuous audible

Loss of Mains supply: Battery powered, audible

Loss of Air or O₂ supply: Pneumatic, audible from blender

System Fail: Visual and intermittent audible

Supplies, Dimensions and Weights

Air and O₂: 3-5 bar

Voltage: 100-120 V 50/60 Hz
220-240 V 50/60 Hz

Power: 120 VA

Fuses: 100-120 V = T500 mA
220-240 V = T200 mA

Protection: Class I Type B

Size, ventilator only: 37 cm W x 31 cm H x 32 cm D

Height on pole: 137 cm

Weight, ventilator only: 10 kg

Designed and manufactured to conform to all relevant international standards for medical devices.

Data subject to change without notice.

tel: +44 (0)20 8681 1414
fax: +44 (0)20 8649 8570
sales@sle.co.uk
www.sle.co.uk

SLE Limited
Twin Bridges Business Park
232 Selsdon Road
South Croydon
Surrey
CR2 6PL
UK



When the smallest thing matters