




DID YOU KNOW – VOLUME SUPPORT

Dräger

VOLUME SUPPORT

What is it?



Volume Support is a spontaneous breathing mode with continuous positive pressure level, including a support in volume.



How is it working?

If the patient is still too weak to manage the complete breathing effort independently, the user can set a target breathing volume and the necessary pressure for it is applied. The set of P_{high} limit alarm avoids too high pressures.

Every inspiratory effort by the patient on the PEEP level that meets the trigger criteria triggers a volume-supported breath. By setting the trigger level, the patient's inspiratory efforts are synchronized.



APN
VTapn520 RRapn 12.0
VS



PC-SIMV

VC-AC

VC-CMV

SPN-CPAP/PS

SPN-CPAP/VS

Other modes

Why is it helpful to improve outcome?

With volume support, the support pressure is automatically adjusted to changes in lung conditions (resistance and compliance) and to the spontaneous breathing demand of the patient.

If the patient's lung mechanics change, the applied pressure varies to keep the set volume constant. The user doesn't have to care about pressure support adjustments, because the device automatically change the pressure to reach the target volume.