

Right. Where you need it. Philips Lumify handheld ultrasound transforms your compatible Android or Apple smart device into a powerful, portable ultrasound system. With advanced processing built directly into the transducer, Lumify delivers exceptional image quality and rapid clinical insights—anytime, anywhere.

See compatible Android devices in the table below and compatible iOS devices in the table on the back.

Lumify Android supported devices*

Device name	Screen size	Connectivity options	Internal storage (GB)	Approximate continuous scan time	USB type	Supports Reacts	Supported OS	Supports "Auto EF" (Yes/No)	Lumify SW version tested
Samsung Galaxy S11 Ultra Tab	14.6 in	Wi-Fi	256	6 h 20 min	C	Yes	16	Yes	5.1.1
Samsung Galaxy A11+ Tab	11 in	Wi-Fi	128	4 h 5 min	C	Yes	16	Yes	5.1.1
Samsung Galaxy Tab Active 5 Pro 5G (SM-X356)	10.1 in	Wi-Fi / 5G	256	5 h 20 min	C	Yes	16	Yes	5.1.1
Samsung Galaxy Active 5 5G Tab (SM-X306)	8 in	Wi-Fi / 5G	128	3 h 25 min	C	Yes	16	Yes	5.1.1
Samsung Galaxy S11 Tablet	11 in	Wi-Fi	128	4 h 16 min	C	Yes	16	Yes	5.1.1
Samsung Galaxy S9 Tab	11 in	Wi-Fi	128	5 h 55 min	C	Yes	16	Yes	5.1.1
Getac ZX10 Gen2 Tablet	10 in	Wi-Fi / LTE	128	6 h 27 min	C	Yes	13	Yes	5.1.1, 3.0
Samsung Galaxy S9 Ultra Tab	14.6 in	Wi-Fi	256	7 h 40 min	C	Yes	16	Yes	5.1.1
Samsung Galaxy S9+ Tab	12.4 in	Wi-Fi	256		C	Yes	16	Yes	5.1.1
Samsung Galaxy S24 Mobile	6.1 in	Wi-Fi / 5G	256	4 h 15 min	C	Yes	16	Yes	5.1.1
Samsung Galaxy S24+	6.7 in	Wi-Fi / 5G	256	6 h	C	Yes	16	Yes	5.1.1
Lenovo Legion Tab 3 (TB321FU)	8.8 in	Wi-Fi	256	4 hr 24 min	C	Yes	15	Yes	5.1
Samsung Galaxy Tab Active 5 Pro (SM-X350) (Wi-Fi Only Tablet)	10.1 in	Wi-Fi	128	6 h 46 min	C	Yes	16	Yes	5.1
Lenovo K10 Gen 2 Tablet (TB311FU)	10.1 in	Wi-Fi	64	2 h 23 min	C	Yes	15	Yes	5.1
Nautix X21 Tablet	6.3 in	Wi-Fi / 5G	64	3 h 35 min	C	Yes	14	Yes	5.1, 3.0
Samsung Galaxy Tab Active 5 (SM-X308)	8 in	Wi-Fi / 5G	128	4 h 15 min	C	Yes	15	Yes	5.1
Samsung Galaxy A9 Tablet	8.7 in	Wi-Fi	64	3 hr 43 min	C	Yes	13	No	4.0.2
Kyocera Tablet (KC-T305CN)	10.1 in	Wi-Fi	64	5 h 31 min	C	Yes	14	No	4.0.2
Getac ZX10	10 in	Wi-Fi / LTE	128	4 hr 39 min	C	Yes	11	No	4.0.1
Samsung Galaxy S23	6.1 in	Wi-Fi / 5G	256	4 h 41 min	C	Yes	14	No	4.0.1
Samsung Galaxy Tab Active 4 Pro	10.1 in	Wi-Fi / 5G	64	4 h 16 min	C	Yes	Out of the box	No	4.0.1
Samsung Galaxy Tab A8	10.5 in	Wi-Fi / 5G	32	3 h 57 min	C	Yes	Out of the box	No	4.0
Huawei Mate 30 Pro 5G	6.53 in	Wi-Fi / 5G	256	3 h 57 min	C	Yes	Out of the box	No	1.9
Huawei Changxiang 2	10.1 in	Wi-Fi	128	3 h 31 min	C	Yes	Out of the box	No	1.7
Panasonic Toughbook (FZ-S1)	7 in	Wi-Fi / LTE	64	2 h 16 min	C	Yes	Out of the box	No	1.9
Panasonic Toughbook (FZ-A3)	10.1 in	Wi-Fi / LTE	128	2 h 19 min	C	Yes	Out of the box	No	1.9

Items in **red** are new to this document revision.

Items in **blue** are compatible, but no longer available for purchase from the manufacturer.

*Devices not included in the preceding table have not undergone compatibility testing with Lumify. Philips cannot guarantee Lumify system performance with non-compatible devices as compatibility testing with all devices on the market is not feasible. The devices have been evaluated using the operating system (OS) versions they were originally equipped with.

Subsequent OS releases, known as maintenance releases, that are higher than the version installed at the time of device launch have not been explicitly verified. Conducting explicit tests for newer major OS releases on older devices, which were initially equipped with a previous OS version, is not feasible. Lumify system functionality and performance may be impacted if significant operating system changes have been applied.

Lumify iOS supported devices*

Device name	Year released	Screen size	Connectivity options	Internal storage (GB)	LPM required	Approx. Continuous Scan Time	Connector type	Supported OS	Lumify SW version tested
iPad Pro 13-inch M5 Chip	2026	13 in	Wi-Fi	256	No	6 h 30 min	C	26.3	5.1.3
iPhone Air	2025	6.5 in	Wi-Fi / 5G	256	Yes	1 h 32 min	C	26.3	5.1.3
iPhone 17	2025	6.27 in	Wi-Fi/5G	256	Yes	1 h 20 min	C	26.2.1	5.1.3
iPad Air 11-inch M3	2025	11 in	Wi-Fi	256	No	NOT TESTED	C	26.2, 18.7.2	5.1.3
iPhone 16	2024	6.1 in	Wi-Fi / 5G	128	Yes	NOT TESTED	C	26.2, 18.7.2	5.1.3
iPad Air 11-inch M2	2024	11 in	Wi-Fi	128	No	NOT TESTED	C	26.2, 18.7.2	5.1.3
iPad Air 13-inch M2	2024	13 in	Wi-Fi	128	No	NOT TESTED	C	26.2	5.1.3
iPad Pro 11-inch (M4)	2024	11 in	Wi-Fi	256	No	NOT TESTED	C	26.3	5.1.3
iPad Pro 13-inch M4	2024	13 in	Wi-Fi	512	No	NOT TESTED	C	26.2, 18.7.2	5.1.3
iPad Pro 12.9-inch 6th Gen	2022	12.9 in	Wi-Fi	128	No	NOT TESTED	C	26.2.1, 18.7.2	5.1.3
iPad Air 5th Gen M1	2022	10.9 in	Wi-Fi	64	No	NOT TESTED	C	26.2, 18.7.2	5.1.3
iPhone 14	2022	6.06 in	Wi-Fi / LTE	128	Yes	NOT TESTED	Lightning	26.1, 18.7.2	5.1.3
iPad 10th Gen	2022	10.9 in	Wi-Fi / LTE	64/256	Yes	NOT TESTED	C	26.3	5.1.3
iPad Mini 6th Gen	2021	8.3 in	Wi-Fi / LTE	64/256	Yes	NOT TESTED	C	26.3	5.1.3
iPhone 17 Pro	2025	6.3 in	Wi-Fi/5G	256	Yes	NOT TESTED	C	26.1	5.1.2
iPhone 16 Pro Max	2024	6.9 in	Wi-Fi/5G	256	Yes	NOT TESTED	C	26.1	5.1.2
iPhone 15	2023	6.12 in	Wi-Fi / LTE	128	Yes	NOT TESTED	C	26.1	5.1.2
iPhone 15 Pro	2023	6.12 in	Wi-Fi / LTE	128	Yes	NOT TESTED	C	26.1, 18.7.2	5.1.2
iPhone 15 Plus	2023	6.69 in	Wi-Fi / LTE	128	Yes	NOT TESTED	C	17.0.3	5.0
iPhone 13	2021	6.1 in	Wi-Fi / LTE	128/25	Yes	NOT TESTED	Lightning	Out of the box	2.0
iPad (9th gen)	2021	10.2 in	Wi-Fi / LTE	64/256	Yes	NOT TESTED	Lightning	Out of the box	2.0

The Items marked in **red** inside the table either newly added or updated to this document revision.

Note: Lumify 2.0 is not compatible with iOS 26.

*Devices not included in the preceding table have not undergone compatibility testing with Lumify. Philips cannot guarantee Lumify system performance with non-compatible devices as compatibility testing with all devices on the market is not feasible. The devices have been evaluated using the operating system (OS) versions they were originally equipped with.

Subsequent OS releases, known as maintenance releases, that are higher than the version installed at the time of device launch has not been explicitly verified. Conducting explicit tests for newer major OS releases on older devices, which were initially equipped with a previous OS version, is not feasible. Lumify system functionality and performance may be impacted if significant operating system changes have been applied.

Please note that Reacts tele-ultrasound solution is not available on Lumify when using the above iOS devices. All iOS devices currently compatible require use of Lumify Power Module and utilize USB-C transducer cables and will continuously scan for 2 to 2.5 hrs.

Devices on the supported list satisfy the following criteria: 1. Passed Philips hardware testing, 2. Is readily available from the device manufacturer.

Philips cannot test every device available on the market, so this list represents a range of devices to satisfy a variety of user needs with respect to form factor, price, battery life and connectivity options.

