



Evita® & More

2021





Mission: "This is what drives us..."

Dear Customer,


Thank you for choosing Dräger. We take pride in providing you with high-quality products and services for all care areas in your hospital. Moreover, our accessories offer you tailor-made solutions for your therapeutic applications. Our original accessories help to protect your patient's safety as ideally as possible. We take just as much care in their specifications, development, and manufacturing, as we do in the manufacturing of our ventilators, anaesthesia workstations, incubators, and monitoring products.

We adhere to strict quality guidelines when it comes to Dräger accessories and consumables. Our high standards of quality have always been at the core of our innovation, and we ensure this by testing all of your products with stringent testing processes in an accredited test center. We don't believe in simply designing and manufacturing technology. Instead, we are driven by supporting you with improved clinical and economical outcomes – that's why we're **your specialist in acute care.**



RESPIRATION PATHWAY

As non-invasive as possible, as invasive as necessary: Discover the variability and diversity of our treatment tools for use along the respiration pathway. They enable administration of protective mechanical ventilation therapy in your ICU.



PREVENT

Avoiding intubation as long as possible helps to decrease the possibility of patients developing ventilator-associated-pneumonia (VAP).

STABILISE

Stabilising your patient is crucial to protecting their lungs from potential conditions such as ventilator-associated lung injury (VALI), which could result in acute lung injury (ALI) or acute respiratory distress syndrome (ARDS).

WEAN

Weaning your patient efficiently reduces incidences of ALI or ARDS.

RECOVER

Ensuring successful weaning and avoiding reintubation improves your patient's recovery.

- A Evita® V800
- B Evita® V600
- C Evita® Infinity® V500*
- D Evita® V300
- E Evita® XL

A

D-5755-2018



B

D-5742-2018



C

D-2106-2015



D

D-2103-2015



E

MT-6127-2008



*Infinity® Acute Care System Workstation Critical Care.

Dräger offers a wide range of software and hardware options for the Evita® family

Compatibility tested accessories and consumables that streamline your workflow and simplify your daily routine are an important part of the Evita® family product portfolio. By providing complete system solutions, Dräger focuses on increasing patient safety and reducing cost of ownership. Please contact your local Dräger sales representative for further information. Not all devices and accessories are available worldwide.

APPLICATIONS FOR RETROFITTING	08	BABYFLOW® PLUS	28
ACCESSORIES FOR BASIC UNIT	10	VENTILATION ACCESSORIES	31
1. Accessories for basic device	10	1. Disposable breathing circuits	31
2. Attachment parts	12	2. Reusable breathing circuits	38
3. Further accessories	12	3. Breathing system filters and HMEs	45
4. Replacement parts for application CO ₂	16	4. Catheter mounts and airway connectors	49
5. Consumables for basic device	17	5. Non-invasive patient interfaces	51
6. Central gas supply	21	6. High-Flow Therapy Nasal Cannula and Accessories	60
BREATHING GAS HUMIDIFIERS	22	7. Pack2Go®	62
1. Humidification with F&P MR 850	22	8. Set2go	63
MEDICATION NEBULIZER	26		
1. Pneumatic medication nebulizer	26		
2. Aerogen Pro	27		

APPLICATIONS FOR RETROFITTING

APPLICATIONS FOR RETROFITTING

Retrofit Kits / Applications

Specify the following when ordering:
Country-specific characteristics

Evita® V800 and V600	8422520
Evita® V300	8420420
IACS Workstation Critical Care Evita® Infinity® V500	8416200

Hardware Applications

Nebulizer

Nurse Call

CO₂ Monitoring

Assembly plate 38 mm (1.5 in) pole
for mounting of ventilation unit

Trolley 2 – 90 cm

Supporting frame with DIN rail

Supporting frame with DIN rail and bail

Cockpit holder

A Gas supply unit GS500

Power supply unit PS500

B Transport supply unit

Transport coupling for Transport supply unit

Gas cylinder holder

C Bed coupling

Bed side rail clamp for bed coupling

Software Applications

D Airway Pressure Release Ventilation

Low Flow PV Loop

AutoFlow/Volume Guarantee

Mandatory Minute Ventilation

Automatic Tube Compensation

Proportional Pressure Support

Monitoring Plus (advanced trend parameter, tabular trends,
graphical trends, data export and advanced loops)

Neonatal Ventilation

E SmartCare®/PS



F Non-Invasive Ventilation

G Smart Pulmonary View

Pressure Support Ventilation

Variable Pressure Support

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



APPLICATIONS FOR RETROFITTING

Evita® XL kits for retrofitting	
Order number for the kits listed below	8414240
CapnoPlus kit	
Nurse call kit	
Evita® Remote kit	
Communication / Evita®Link kit	
NeoFlow kit	
Mask ventilation (NIV) kit	
Evita® XL 5.10 Evita® XL software kit	
6.0 Evita® XL software kit	
A SmartCare® kit; not for USA	
SmartCare® Update 1.1	
SmartCare® Capno package kit	
Lung protection package kit	
O ₂ -therapy kit	
Proportional Pressure Support (PPS) kit; not for USA	
Software 7.06, NIV Plus	
Plug for nurse call connection	1846248
Accessories for nurse call option.	
B CO₂ mainstream sensor	6871500
Replacement parts for SmartCare® option (not for USA).	

For patient application see O2Star Portfolio.

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



MT-0619-2008

B



MT-0986-2007

ACCESSORIES FOR BASIC UNIT

ACCESSORIES FOR BASIC UNIT

1. ACCESSORIES FOR BASIC DEVICE

ACCESSORIES FOR TROLRIES

EVITA® INFINITY® V500 AND EVITA® V300

A	Gas cylinder holder	G93110				
B	Bedside rail clamp for bed coupling	8419801				
C	Rail clamp to mount F&P* humidifier	8411074				

ACCESSORIES FOR MEDICAL COCKPIT NOT MOUNTED ON TROLLEY

Module M26 for mounting of Infinity® Medical Cockpit
Specify the following when ordering:
Country-specific characteristics

		G96500				
	System cable, length 0.70 m (27.6 in)	MS18782				
	System cable, length 2.80 m (9.2 ft)	MS18781				

ACCESSORIES FOR TROLRIES

D	Pressure reducer OxyLine, customized configuration of pressure reducers	MP00100				
	OxyLine's modular design means you can accurately choose the pressure reducer you need from the following components:					
	– O₂ or AIR pressure reducer with 200 bar – Manometer, vertical or horizontal – Flow regulator settings (1, 2, 3, 4, 6, 8, 12 and 15 L/min) and 9/16" sterile water connector – Flowmeter for up to 16 L/min with connector for Dräger reusable humidifier, sterile water connector and 6 mm (0.2 in) nozzle for an oxygen mask – Self-closing NIST connector – Centralized gas supply coupling, e.g. to connect up a mains-operated ventilation device – Ejector for short-time suction of eg. secret					
E	O₂-NIST/NIST Connecting hose 1.2 m (3.9 ft)	M36053				
	A choice of additional hoses for connecting the pressure reducer with the Evita® XL: 8601697.					

*Produced by Fisher & Paykel

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

TROLLEY AND ACCESSORIES FOR TROLLEY EVITA® XL

A	EvitaMobil with extension (USA)	8414455
----------	--	----------------

ACCESSORIES FOR TROLRIES EVITA® XL

B	Bottle holder set, EvitaMobil	8411970
----------	--------------------------------------	----------------

C	DC battery 12 V/17 Ah (box of 1; 2 required)	1843303
----------	---	----------------

	Battery cable kit for EvitaMobil	8411822
--	---	----------------

COMPRESSOR EVITA® XL

D	Breathing air compressor for supply with compressed air, can be integrated in the EvitaMobil	8414350
----------	---	----------------

	Standalone base for compressor	8411520
--	---------------------------------------	----------------

	Air hose, 0.8 m, NIST neutral	8411517
--	--------------------------------------	----------------

	Air hose, 0.8 m, NIST ISO 32	8411538
--	-------------------------------------	----------------

	Air hose, 5 m, NIST ISO 32	8411541
--	-----------------------------------	----------------

	Air hose, 5 m, NIST neutral	8411519
--	------------------------------------	----------------

	Air hose, 0.8 m, DISS	8414454
--	------------------------------	----------------

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

2. ATTACHMENT PARTS

A	Instrument tray, complete (EvitaXL only)	8414828
	Side holder set	8411615

STORAGE

B	Compact rail, for attaching the display unit to the 38 mm (1.5 in) pole	2M85337
C	CS connection hose holder	2M85446
D	Hook	M26349
	Hook for accessories, e.g. breathing bag.	

3. FURTHER ACCESSORIES

	Coarse dust filter for ventilator, for patient part heating, as of 2005	8415978
	Flow sensor flap; for patient part heating, as of 2005	8415973
	ILV cable, Evita® 4/XL	8411794
	ILV cable, Evita® 1/2	8411793
	Monitor bracket kit, Evita®, with trolley	8415787
	Monitor bracket kit, Evita®, without trolley	8415850
	Printer cable	8306489

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

MEDIBUS cable	8416326
Software Protocol, MEDIBUS for V and VN ventilators	9039527
Software Protocol, MEDIBUS.X for V and VN ventilators	9052608
USB flash drive	8416347
NO THERAPY	
NO application device*	8421300
Specify the following when ordering:	
Country-specific characteristics	

*Produced by EKU

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

CONTENT

RETROFITTING

BASIC UNIT

HUMIDIFIERS

MEDICATION NEBULIZER

BABYFLOW® PLUS

VENTILATION

ACCESSORIES FOR BASIC UNIT

A	Test lung SelfTestLung, max. 1 L (33.81 fl. oz.), reusable, not made with natural rubber latex	MP02400
	Spare parts set SelfTestLung, with locking brace and attachment bolt	MP02401
B	Test lung NeoTestLung, disposable, N, bellows, 10 pcs., not made with natural rubber latex	MP02425
C	Demo lung to simulate neonatal lungs, with suitcase	8418077

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

HINGED ARMS

Keeping everything close at hand. To leave your hands free to care for the patient.

- | | |
|--|-----------------------|
| <p>A Hinged arm Infinity® ACS, for positioning and fixation of breathing hoses, this arm can be fixed into place without a hand movement through a patent registered friction torque-breaking system, length 1,080 mm, max. load 0.5 kg, includes hose clamps: MP00693 (1×), 8417396 (2×)</p> | <p>MP00690</p> |
| <p>B Hinged arm 240°, 3 hinges, fixable, ball holder system, length 1,155 mm, max. load 0.3 kg, includes hose clamps: MP00493 (1×), 8409841 (1×), 8409746 (1×)</p> | <p>8409609</p> |
| <p>C Hinged arm Quickstop 2, with universal ball holder system for optional accessories length 858 mm, max. load 0.4 kg, includes hose clamps: MP00493 (1×), 8404761 (2×), 8409841 (1×)</p> | <p>2M85706</p> |

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

4. REPLACEMENT PARTS FOR APPLICATION CO₂

A	CO ₂ mainstream sensor	6871950
B	CapnoSmart CO ₂ main flow sensor Replacement parts for SmartCare® option (not for USA).	6871500
CO ₂ CUVETTES, DISPOSABLE		
C	① CO ₂ cuvette disposable, pediatric, 10 pcs.	MP01063
	② CO ₂ cuvette disposable, adult, 10 pcs.	MP01062
CO ₂ CUVETTES, REUSABLE		
D	① CO ₂ cuvette reusable, pediatric	6870280
	② CO ₂ cuvette reusable, adult	6870279
ACCESSORIES		
	Calibration set, with test connection and test gas cylinder	8412710
	Test gas cylinder 5 Vol.% CO ₂ , 95 Vol. % N ₂	6850435
E	Parking holder for CO ₂ sensor	8416352
	Parking holder for CO ₂ sensor	8412840

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

5. CONSUMABLES FOR BASIC DEVICE

A	Expiratory valve for Evita® 2 dura, Evita® 4 edition, Evita® XL, disposable, 10 pcs.	8414776
	Expiratory valve for Evita® 2 dura, Evita® 4 edition, Evita® XL, reusable	8410580
B	Expiratory valve for Evita® 2 dura, Evita® 4 edition, Evita® XL, disposable, 10 pcs.	MP02600
C	Infinity® ID expiratory valve, disposable, 10 pcs.	MP01060
D	Expiratory valve for Evita® Infinity® V500, reusable	8416750
E	Expiratory valve for Evita® V800, V600 and Savina® 300, reusable	8417050
F	Expiratory valve, disposable, 10 pcs.	MP01061
	Expiratory valve, for neonatal/pediatric patients (Non-RFID), reusable	8415296
G	Infinity® ID expiratory valve, for neonates/pediatric patients, reusable	8415270

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



D-54701-2012

B



D-54687-2012

C



D-5887-2018

D



D-32547-2011

E



D-45109-2015

F



D-11898-2019

G



D-45107-2015

ACCESSORIES FOR BASIC UNIT

	Membrane for Evita® 2 dura, Evita® 4, Evita® XL, for expiratory valve	8412015
A	Water trap for Evita® 2 dura, Evita® 4, Evita® XL, reusable	8413125
B	Infinity® ID Expiratory filter, disposable, 20 pcs.	MP01780
C	Expiratory filter, disposable, 20 pcs.	MP01781
	Inspiratory valve, reusable	8417015
	Muffler	8415230
COMPATIBLE WITH THE EXPIRATORY VALVE FOR ADULTS/PEDIATRIC PATIENTS		
	Replacement diaphragm for reusable expiratory valve	8413661
	Infinity® ID Inspiratory valve, reusable	8415195
	Ambient air filter/Filter for fan, set of 10 pcs.	8416349
	Ambient air filter/Filter for fan, set of 10 pcs.	8412384

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



MT-220-99

B



D-12244-2010

C




D-43830-2012

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

Hot-wire anemometry: Accurate and intelligent flow measurement

Spirolog®, neonatal flow sensors and SpiroLife® are based on proven hot-wire anemometry techniques. Due to extremely fast and accurate measurement results, as well as extremely low inspiratory and expiratory resistances, it is clearly superior to other methods. You only need one flow sensor for all patients. Neonatal flow sensors offer customized solutions for the sole exception of the specific needs of neonates in intensive care. The sterilizable neonatal flow sensors and SpiroLife® combine the advantages of hot-wire technology with the possibility of sterilizability. Therefore, they integrate seamlessly into your hospital's sterilization processes.

A	<u>Flow sensor Spirolog®, PC/ABS, 5 pcs.</u>	<u>8403735</u>
B	<u>Flow sensor SpiroLife®, PPSU, sterilizable</u>	<u>MK01900</u>
C	<u>Flow sensor Infinity® ID, PC/ABS, 5 pcs.</u>	<u>6871980</u>
D	<u>Neonatal flow sensor with Y-piece, PSU, sterilizable</u>	<u>8410185</u>
E	<u>Neonatal flow sensor ISO 15, PSU, sterilizable</u>	<u>8411130</u>
F	<u>Insert for neonatal flow sensor, PSU, sterilizable, 5 pcs.</u>	<u>8410179</u>
	Flow sensor cover	8417163



ACCESSORIES FOR BASIC UNIT

OXYGEN SENSOR BY DRÄGER

Oxygen for Life

Oxygen sensors have the important task of measuring and monitoring the inspiratory and/or expiratory O₂ concentration in anesthesia and respiratory devices. Oxygen sensors from Dräger satisfy the highest national and international requirements. Furthermore, they have been extensively tested and fine-tuned for proper functional performance with Dräger devices and systems. Dräger's years of experience in the field of sensors, and their high quality has even convinced NASA¹⁾. Dräger O₂ sensors are electrochemical transducers used to measure O₂ partial pressure. They operate on the principle of a galvanic cell.

Examples of locations where the sensors are used include the Fabius® family, Evita® family and the Primus®.

To take into account the high value of oxygen measurement, we only recommend using O₂ sensors from Dräger.

A	O₂ sensor capsule, for mainstream measurement	6850645
----------	---	----------------

CABLES

B	Flow sensor cable, for neonatal flow sensor	8409626
C	Flow sensor cable for V500 with Neo option and VN500	8416600

¹⁾Dräger, Die Geschichte des Unternehmens, 9070230, page 28

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ACCESSORIES FOR BASIC UNIT

6. CENTRAL GAS SUPPLY

A	Dräger supply hoses, selectable by the adjacent configuration number	8601697
	Design, color coding and connection thread of all Dräger supply hoses comply with the respective legal regulations in the indicated countries.	

FILTER FOR CENTRAL SUPPLY CONNECTING HOSES

High Pressure Line Filter for medical air	MP03901
High Pressure Line Filter for oxygen	MP03902
Filter element with O-ring	MP03903
Filter bowl with O-ring	MP03904

ADAPTER HOSES

Adapter O₂ NIST for Australia	M34446
Adapter AIR NIST for Australia	M34448

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BREATHING GAS HUMIDIFIERS

BREATHING GAS HUMIDIFIERS

1. HUMIDIFICATION WITH F&P* MR 850

Humidifier basic unit F&P MR 850

Specify the following when ordering:

Country-specific characteristics

8414144

- A Temperature probe/flow sensor for F&P MR 850,
for all breathing circuits, 1.1 m (3.6 ft)

8414989

- B Temperature probe/flow sensor for F&P MR 850,
for all breathing circuits, 1.5 m (4.92 ft)

8414966

- C Hose heater adapter for F&P MR 850,
for reusable breathing circuit

8414968

- D Hose heater adapter for F&P MR 850,
for disposable breathing circuit

8414992

- E Hose heater adapter for F&P MR 850,
for double-heated disposable breathing circuit

MP01006

- F Draw wire for hose heater 1.5 m,
for reusable breathing circuit

8411050

*Produced by Fisher & Paykel

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-1766-2019



B

D-1767-2019



C

D-1766-2019



D

D-1768-2019




E

MT-4048-2006



F

MT-142-98



BREATHING GAS HUMIDIFIERS

REQUIRED ACCESSORIES FOR F&P MR 850

A	Humidifier chamber F&P* MR 290, disposable, 10 pcs.	8418282
B	Humidifier chamber F&P* MR 340, reusable, N/P	8411047
C	Humidifier chamber F&P* MR 370, reusable, A	8412217
D	Humidifier chamber VentStar®, disposable, 10 pcs.	MP00234
E	Humidifier chamber VentStar®, disposable, N, 10 pcs.	MP02590
F	Humidifier chamber, 10 pcs.	MP02605

*Produced by Fisher & Paykel

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A	
B	
C	
D	
E	
F	

MT-0099-2004

D-7591-2019

D-7593-2019

D-24199-2009

D-42729-2015

D-94685-2013


BREATHING GAS HUMIDIFIERS

HOLDER

A	Rail clamp to mount F&P humidifier*	8411074
B	Swivel-mounted humidifier holder for trolley	G93111
C	Holder for humidifier, lateral standard rail	8416325
D	Universal holder with standard rail	G93140
E	Holder for F&P humidifier U-shaped, for trolley, with rail clamp*	MP01002
	Humidifier holder	8411956
	Holder, lateral (F&P 900 MR 309), including 8411074*	MP01009
F	Water bag pole (F&P 800 MR 209), holder for water bag, length 70 cm (27,6 in)*	MP01008

*Produced by Fisher & Paykel

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BREATHING GAS HUMIDIFIERS

HOSES*

A	Hose heater, A, 1.1 m (43,3 in)	8411045
B	Hose heater for F&P MR 700/730 and MR 850, for reusable breathing circuit, 1.1 m (43,3 in)	8414993

*Produced by Fisher & Paykel

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-2264-2004



B

MT-2266-2004



MEDICATION NEBULIZER

MEDICATION NEBULIZER
1. PNEUMATIC MEDICATION NEBULIZER

A	Medication nebulizer for Evita® 2 dura, Evita® 4 edition, Evita® XL	8412935
B	Connector set for medication nebulizer, for neonatal breathing circuit	8411031

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



MEDICATION NEBULIZER

2. AEROGEN PRO*

Aerogen Pro works with an electronic micropump. The nebulizer can be used with Evita® 2 dura, Evita® 4 edition and Evita® XL.

Aerogen Pro is ideally suited to nebulizer therapy requirements for inline mechanical ventilation and it accommodates all patient age ranges from neonates to adults. It operates without changing patient ventilator parameters, without heating or degrading medications and it can be refilled without interrupting ventilation.

A	Medication nebulizer Aerogen Pro	MP01010						
B	Nebulizer unit with filler cap for Aerogen Pro	MP01030						
C	T-piece with plug for Aerogen Pro, A, ISO 22 mm (0.87 in), 5 pcs.	MP01031						
D	T-piece set with plug for Aerogen Pro, N, 15 mm (0.59 in), w/ adapter 10 mm (0.4 in), 5 pcs.	MP01032						
E	T-connector for Aerogen Pro, 12 mm (0.47 in), with plug, 5 pcs.	MP01037						
F	Controller cable for Aerogen Pro	MP01036						
G	Filler cap for Aerogen Pro, silicone, 5 pcs.	MP01035						

*Produced by Aerogen

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BABYFLOW® PLUS

BABYFLOW® PLUS

BabyFlow® plus – The superior non-invasive respiratory support system, designed to maximise performance and patient comfort.

BABYFLOW® PLUS

A	<u>BabyFlow® plus, disposable, size S, 20 pcs.</u>	<u>MP03700</u>				
	<u>BabyFlow® plus, disposable, size M, 20 pcs.</u>	<u>MP03701</u>				
	<u>BabyFlow® plus, disposable, size L, 20 pcs.</u>	<u>MP03702</u>				
	<u>BabyFlow®, disposable, 20 pcs.</u>	<u>8418583</u>				
B	<u>Demo package BabyFlow® plus, disposable, size S</u>	<u>MP03703</u>				
	<u>Demo package BabyFlow® plus, disposable, size M</u>	<u>MP03704</u>				
	<u>Demo package BabyFlow® plus, disposable, size L</u>	<u>MP03705</u>				
C	<u>Soft connector N, disposable, 20 pcs.</u>	<u>MP03826</u>				

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-7477-2016



B

D-2195-2015



C

D-42601-2015



BABYFLOW® PLUS

STANDARD PRONGS

Soft anatomically formed nasale silicone prongs, disposable.

A	① Prong BabyFlow®, disposable, size XS, 10 pcs.	8418415						
	② Prong BabyFlow®, disposable, size S, 10 pcs.	8418605						
	③ Prong BabyFlow®, disposable, size M, 10 pcs.	8418416						
	④ Prong BabyFlow®, disposable, size L, 10 pcs.	8418531						
	⑤ Prong BabyFlow®, disposable, size XL, 10 pcs.	8418417						
B	Prong BabyFlow®, disposable, size XXL, 10 pcs.	8418617						

WIDE PRONGS



	Prong BabyFlow®, wide, disposable, size 3, 10 pcs.	8418603						
	Prong BabyFlow®, wide, disposable, size 4, 10 pcs.	8418607						
	Prong BabyFlow®, wide, disposable, size 5, 10 pcs.	8418609						
	Prong BabyFlow®, wide, disposable, size 6, 10 pcs.	8418611						

MASKS

Nasal masks of silicone and TPE complete the portfolio and increase the comfort for the patient.

C	① Neo-mask, silicone, disposable, size S, 10 pcs.	8418491						
	② Neo-mask, silicone, disposable, size M, 10 pcs.	8418490						
	③ Neo-mask, silicone, disposable, size L, 10 pcs.	8418619						
D	BabyFlow® Mask, size 4, Box 10	MP01494						
E	BabyFlow® Mask, size 5, Box 10	MP01495						
F	BabyFlow® Mask, size 6, Box 10	MP01496						

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



BABYFLOW® PLUS

HEADGEARS

The new headgear enable the customer to choose the most comfortable fit and combination for each of their tiny patients.

A	Headgear, disposable, size XXS, 17 – 20 cm (6.7 – 7.9 in) head circumference, 5 pcs.	MP03880
	Headgear, disposable, size XS, 20 – 23 cm (7.9 – 9.1 in) head circumference, 5 pcs.	MP03881
B	Headgear, disposable, size S, 23 – 27 cm (9.1 – 10.6 in) head circumference, 5 pcs.	MP03882
	Headgear, disposable, size M, 27 – 31 cm (10.6 – 12.2 in) head circumference, 5 pcs.	MP03883
	Headgear, disposable, size L, 31 – 37 cm (12.2 – 14.6 in) head circumference, 5 pcs.	MP03884
	Headgear, disposable, size XL, 37 – 42 cm (14.6 – 16.5 in) head circumference, 5 pcs.	MP03885

CAPS


The disposable and breathable microfiber caps (95% cotton, 5% spandex) available in seven sizes (color-coded).

C	① Cap BabyFlow®, disposable, size XS, white, 18 – 20.5 cm (7.1 – 8.1 in) head circumference, 5 pcs.	8418533
	Cap BabyFlow®, disposable, size S, yellow, 20.5 – 23 cm (8.1 – 9.1 in) head circumference, 5 pcs.	8418534
	Cap BabyFlow®, disposable, size M, red, 23 – 25.5 cm (9.1 – 10 in) head circumference, 5 pcs.	8418535
	② Cap BabyFlow®, disposable, size L, green, 25.5 – 28 cm (10 – 11 in) head circumference, 5 pcs.	8418536
	③ Cap BabyFlow®, disposable, size XL, orange, 28 – 30.5 cm (11 – 12 in) head circumference, 5 pcs.	8418537
	Cap BabyFlow®, disposable, size XXL, blue, 30.5 – 33 cm (12 – 13 in) head circumference, 5 pcs.	8418538
D	Cap BabyFlow®, disposable, size XXL+, white/black, 33 – 36 cm (13 – 14.2 in) head circumference, 5 pcs.	8418539

FIXATION STRAPS

Fixation strap BabyFlow®, disp., 13 cm (5,12 in), 10 pcs.	MP03870
Fixation strap BabyFlow®, disp., 15 cm (5.91 in), 10 pcs.	MP03871
Fixation strap BabyFlow®, disp., 18 cm (7,09 in), 10 pcs.	MP03872
Fixation strap BabyFlow®, disp., 20 cm (7,87 in), 10 pcs.	MP03873

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



VENTILATION ACCESSORIES

VENTILATION ACCESSORIES

1. DISPOSABLE BREATHING CIRCUITS

INFINITY® ID PRODUCTS FOR ADULTS

A	Breathing circuit Infinity® ID, disposable, basic, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP01350
B	Breathing circuit Infinity® ID, disposable, basic, 2.5 m (8.2 ft), 10 pcs., not made with natural rubber latex	MP01348
C	Breathing circuit Infinity® ID, disposable, water trap, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP01326
D	Breathing circuit Infinity® ID, disposable, coaxial, 1.6 m (5.25 ft), 20 pcs., not made with natural rubber latex	MP00342
E	Breathing circuit Infinity® ID, disposable, smoothbore, 1.6 m (5.25 ft), 25 pcs., not made with natural rubber latex	MP00340
F	Breathing circuit Infinity® ID, disposable, smoothbore, 3 m (9.8 ft), 10 pcs., not made with natural rubber latex	MP01325
G	Breathing circuit Infinity® ID, disposable, extendable, 2.2 m (7.22 ft), 25 pcs., not made with natural rubber latex	MP01347

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



DISPOSABLE BREATHING CIRCUITS

A Breathing circuit Infinity® ID,
disposable, coaxial, 1.9 m (6.2 ft), 20 pcs.,
not made with natural rubber latex MP03375

B Breathing circuit Infinity® ID,
disposable, coaxial, 2.3 m (7.55 ft), 20 pcs.,
not made with natural rubber latex MP03377

C Breathing circuit Infinity® ID,
disposable, coaxial, with water traps,
2.1 m (6.89 ft), 10 pcs.,
not made with natural rubber latex MP00378

INFINITY® ID PRODUCTS FOR PEDIATRICS

D Breathing circuit Infinity® ID,
disposable, P, basic, 1.8 m (5.91 ft), 25 pcs.,
not made with natural rubber latex MP01340

E Breathing circuit Infinity® ID,
disposable, P, water trap, 1.8 m (5.91 ft), 25 pcs.,
not made with natural rubber latex MP01341

INFINITY® ID PRODUCTS FOR NEONATES

F Breathing circuit Infinity® ID,
disposable, N, heated, 1.2 m (3.9 ft), 10 pcs.,
not made with natural rubber latex MP00343

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL




DISPOSABLE BREATHING CIRCUITS

PRODUCTS WITHOUT INFINITY® ID FOR ADULTS

A	Breathing circuit Vent Set, disposable, basic, 1.5 m (4.92 ft), 25 pcs., not made with natural rubber latex	MP00300				
B	Breathing circuit VentStar®, disposable, basic, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP00349				
C	Breathing circuit VentStar®, disposable, basic, 2.5 m (8.2 ft), 10 pcs., not made with natural rubber latex	MP00350				
D	Breathing circuit Vent Set, disposable, extendable, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP00305				
E	Breathing circuit VentStar®, disposable, extendable, 2.2 m (7.22 ft), 25 pcs., not made with natural rubber latex	MP00355				

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



DISPOSABLE BREATHING CIRCUITS

A	Breathing circuit Vent Set, disposable, coaxial, 1.5 m (4.92 ft), 20 pcs., not made with natural rubber latex	MP00315
B	Breathing circuit Vent Set, disposable, coaxial, with holder, 1.5 m (4.92 ft) 20 pcs., not made with natural rubber latex	MP00316
C	Breathing circuit VentStar®, disposable, water trap, w/o LL, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP00334
D	Breathing circuit VentStar®, disposable, water trap, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP00337
E	Breathing circuit VentStar®, disposable, smoothbore, for MRI, 3 m (9.8 ft), 10 pcs., not made with natural rubber latex	MP00327
	VentStar® MRI, 3 m (9,8 ft), 10 pcs.	MP04800
	Breathing circuit VentStar®, disposable, coaxial, 1.7 m, 20 pcs., not made with natural rubber latex	MP03373
F	Breathing circuit VentStar®, disposable, coaxial, 1.9 m (6.2 ft), 20 pcs., not made with natural rubber latex	MP03374
	Breathing circuit VentStar®, disposable, coaxial, 2.3 m, 20 pcs., not made with natural rubber latex	MP03376
G	Breathing circuit VentStar®, disposable, coaxial, with water traps, 2.1 m (6.89 ft), 10 pcs., not made with natural rubber latex	MP00373

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL




DISPOSABLE BREATHING CIRCUITS

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A	VentStar® Helix, dual heated, 10 pcs.	MP02606
B	VentStar® Helix, heated, 10 pcs.	MP02607
C	Humidifier chamber, 10 pcs.	MP02605
	Breathing circuit Fisher & Paykel RT212, disposable, inspiratory heated, 10 pcs.*	8414991

*Produced by Fisher & Paykel




DISPOSABLE BREATHING CIRCUITS

DISPOSABLE BREATHING CIRCUITS FOR PEDIATRICS

A	Breathing circuit VentStar®, disposable, P, basic, w/o LL, 1.5 m (4.92 ft), 25 pcs., not made with natural rubber latex	MP00338
B	Breathing circuit VentStar®, disposable, P, basic, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP00351
C	Breathing circuit VentStar®, disposable, P, water trap, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP00361
D	Breathing circuit VentStar®, disposable, P, water trap, w/o LL, 1.8 m (5.91 ft), 25 pcs., not made with natural rubber latex	MP00362

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



DISPOSABLE BREATHING CIRCUITS

DISPOSABLE BREATHING CIRCUITS
FOR NEONATES

A	Breathing circuit VentStar®, disposable, N, heated, 1.2 m (3.9 ft), 10 pcs., not made with natural rubber latex	MP00308
B	Breathing circuit VentStar®, disposable, N, heated, w/o chamber, 1.2 m (3.9 ft), 10 pcs., not made with natural rubber latex	MP00309
C	Breathing circuit VentStar® Helix, disposable, N, dual heated, 1.7 m (5.58 ft), 10 pcs., not made with natural rubber latex	MP02650
D	Breathing circuit VentStar® Helix, disposable, N, heated, 1.7 m (5.58 ft), 10 pcs., not made with natural rubber latex	MP02608
E	Double connector 11 ID/15 OD, is required when using MP02650 and MP02608, set of 20 pcs.	MP02655
F	Humidifier chamber VentStar®, disposable, 10 pcs.	MP00234
G	Humidifier chamber VentStar®, disposable, N, 10 pcs.	MP02590
	Breathing circuit Fisher & Paykel RT225, disposable, N, inspiratory heated, 10 pcs.*	8414994

*Produced by Fisher & Paykel

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL


REUSABLE BREATHING CIRCUITS

2. REUSABLE BREATHING CIRCUITS

BLUESET VENTILATION

A	BlueSet Ventilation HME (P/A)	8412860
	Consists of:	
1	Breathing hose, silicone, reusable, A, 1.5 m (4.92 ft), connector 22 mm/22 mm (0.87/0.87 in)	2166046
2	Y-piece, reusable, A	M25650


EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



REUSABLE BREATHING CIRCUITS

Silicone breathing hose (adult) with standard sleeves 22 mm/22 mm (0.9 in/0.9 in)	
A Breathing hose, silicone, reusable, 0.35 m (13.7 in), connector 22 mm/22 mm (0.87/0.87 in)	2165619
B Breathing hose, silicone, reusable, 60 cm (1.97 ft), connector 22 mm/22 mm (0.87/0.87 in)	2165627
C Breathing hose, silicone, reusable, 1.1 m (3.6 ft), connector 22 mm/22 mm (0.87/0.87 in)	2166038
D Breathing hose, silicone, reusable, 1.5 m (4.9 ft), connector 22 mm/22 mm (0.87/0.87 in)	2166046

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL





REUSABLE BREATHING CIRCUITS

Silicone breathing hose (pediatric) with standard sleeves with one 10 mm (0.4 in) and one 22 mm (0.9 in) sleeve and 10 mm (0.4 in) inner diameter

A	Breathing hose, silicone, reusable, P, 40 cm (15.8 in), connector 10 mm/22 mm, Ø 10 mm	2165856
B	Breathing hose, silicone, reusable, P, 60 cm (1.97 ft), connector 10 mm/22 mm, Ø 10 mm	2165821
C	Breathing hose, silicone, reusable, P, 1.1 m (3.61 ft), connector 10 mm/22 mm, Ø 10 mm	2165651
D	Breathing hose, silicone, reusable, P, 1.5 m (4.92 ft), connector 10 mm/22 mm, Ø 10 mm	2165961

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



REUSABLE BREATHING CIRCUITS

Silicone breathing hose (pediatric) with standard sleeves 10 mm (0.4 in) on both sides and 10 mm (0.4 in) inner diameter

A Breathing hose, silicone, reusable, P, 60 cm (1.97 ft), connector 10 mm/10 mm, Ø 10 mm

2165848


Silicone breathing hose (pediatric) with standard sleeves 22 mm (0.9 in) on both sides and 10 mm (0.4 in) inner diameter

B Breathing hose, silicone, reusable, P, 1.1 m (3.61 ft), connector 22 mm/22 mm, Ø 10 mm

2165678

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



0-36215-2011

B



0-36213-2011

VENTILATION ACCESSORIES

Y-PIECES AND CONNECTORS

A	Y-piece, reusable, for temperature probe adapter 900 MR 532	8412888
B	Temperature probe adapter 900 MR 532, for F&P* humidifiers	MP01007
C	Y-piece, reusable, for Dräger temperature probe	8405435
D	Y-piece, reusable, P, angled, with suction opening	8403075
	Cap, reusable, 7 mm (0.28 in)	8401645
E	Elbow, reusable	M25649

*Produced by Fisher & Paykel

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-1321-2008



B

MT-4036-2006



C

MT-2688-2004



D

MT-2728-2004



E

MT-2688-2004



VENTILATION ACCESSORIES

A	ISO double connector, reusable, 22 mm/22 mm (0.87/0.87 in)	M25647
	Connector for breathing hose, reusable, A	8403076
	T-piece connector, reusable	8411044
	Double connector, outer diameter 11 mm (0.43 in)	8409897
B	Water trap, for use between the ventilator hoses	8404985
C	Water trap, reusable, (N/P)	8409627

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



VENTILATION ACCESSORIES

	<u>Set of caps, set of 5 pcs.</u>	<u>8402918</u>
	<u>Corrugated hose</u>	<u>8402041</u>
	<u>Corrugated hose, flexible, length 0.14 m (5.5 in)</u>	<u>8403333</u>
A	<u>Corrugated hose, flexible, length 0.13 m (5.3 in)</u>	<u>8409634</u>
	<u>Hose clip</u>	<u>8403566</u>
B	<u>ISO double connector, reusable, 22 mm/22 mm (0.87/0.87 in)</u>	<u>M25647</u>
C	<u>Temperature sensor</u>	<u>8405371</u>

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-1317-2020



B

MT-0732-2008



C

D-2368-2011



BREATHING SYSTEM FILTERS AND HMEs

3. BREATHING SYSTEM FILTERS AND HMEs

ELECTROSTATIC FILTERS AND HMEs

The TwinStar® filter/HME from Dräger is a combination of an efficient HME and a high-performance breathing system filter, thereby significantly contributing to infection prophylaxis in ventilation treatment.

- Combination of efficient HME and high-performance breathing system filter
- High retention rates for bacteria and viruses
- High return of heat and moisture
- Reduced condensation in the breathing circuit
- Minimal workload and low cost
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and blue color coding for quick identification

A	Electrostatic filter and HME TwinStar® 90, disposable, 90 mL (3.04 fl. oz.) deadspace, 50 pcs.	MP01800
B	Electrostatic filter and HME TwinStar® 65A, disposable, 65 mL (2.2 fl. oz.) deadspace, angled, 50 pcs.	MP01810
C	Electrostatic filter and HME TwinStar® 55, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01805
D	Electrostatic filter and HME TwinStar® 25, disposable, 25 mL (0.85 fl. oz.) deadspace, 50 pcs.	MP01815
E	Electrostatic filter and HME TwinStar® 10A, disposable, 10 mL (0.34 fl. oz.) deadspace, angled, 50 pcs.	MP01825
F	Electrostatic filter and HME TwinStar® 8, disposable, 8 mL (0.27 fl. oz.) deadspace, 50 pcs.	MP01820
G	Mechanical filter and HME TwinStar® HEPA, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01801

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

MT-4385-2007	A	
MT-4745-2007	B	
MT-1164-2006	C	
MT-0105-2008	D	
MT-4397-2007	E	
MT-4396-2007	F	
MT-4388-2007	G	

BREATHING SYSTEM FILTERS AND HMEs

ELECTROSTATIC FILTERS

The CareStar® breathing system filter significantly contributes to infection prophylaxis in anesthesia and ventilation treatment through high retention rates.

- High-performance electrostatic filter
- Patient-side or device-side use
- Minimal workload and low cost
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and red color coding for quick identification

A	Electrostatic filter CareStar® 45, disposable, 45 mL (1.52 fl. oz.) deadspace, 50 pcs.	MP01755					
B	Electrostatic filter CareStar® 40A, disposable, 40 mL (1.35 fl. oz.) deadspace, angled, 50 pcs.	MP01765					
C	Electrostatic filter CareStar® 30, disposable, 30 mL (1.01 fl. oz.) deadspace, 50 pcs.	MP01770					

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-4387-2007



B

MT-4747-2007



C

MT-0011-2007



BREATHING SYSTEM FILTERS AND HMEs

MECHANICAL FILTERS

The SafeStar® mechanical high performance breathing system filter shows excellent retention rates, thereby significantly contributing to infection prophylaxis in anesthesia and ventilation treatment.

- High-performance mechanical HEPA filter
- High hydrophobicity
- Patient-side or device-side use
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and red color coding for quick identification

A	Mechanical filter SafeStar® 80, disposable, 80 mL (2.71 fl. oz.) deadspace, 50 pcs.	MP01785					
B	Mechanical filter SafeStar® 60A, disposable, 60 mL (2.03 fl. oz.) deadspace, angled, 50 pcs.	MP01795					
C	Mechanical filter SafeStar® 55, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01790					

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

MT-4386-2007



B

MT-4746-2007



C

MT-1165-2006



BREATHING SYSTEM FILTERS AND HMEs

HMEs

HumidStar® heat and moisture exchangers from Dräger partially take over the function of the upper airways and efficiently heat and humidify the inspired air.

- High return of heat and moisture
- Reduced condensation in the breathing circuit
- Minimal workload and low cost
- Sampling connector with tethered cap
- Sampling port in convenient 45° angle
- Transparent housing for visual control
- Standardized connectors for safe connection to other components
- Clear labeling and green color coding for quick identification

A	HME HumidStar® 55, disposable, 55 mL (1.86 fl. oz.) deadspace, 50 pcs.	MP01730
B	HME HumidStar® 25, disposable, 25 mL (0.85 fl. oz.) deadspace, 50 pcs.	MP01735
C	HME HumidStar® 10A, disposable, 10 mL (0.34 fl. oz.) deadspace, angled, 50 pcs.	MP01740
D	HME HumidStar® 2 Plus, disposable, 2 mL (0.07 fl. oz.) deadspace, 100 pcs.	MP05845
E	HME HumidStar® Trach plus, disposable, 6 mL deadspace, 100 pcs.	MP05750

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



CATHETER MOUNTS AND AIRWAY CONNECTORS

4. CATHETER MOUNTS AND AIRWAY CONNECTORS

- Smoothbore products provide low resistance to flow while external coil reduces kinks and blockages
- Extendable product provides a lightweight and flexible solution
- Double swivel connector increases ease of use
- Double cap to maintain PEEP during suctioning or bronchoscopy
- Each catheter mount and airway connector is equipped with a red safety cap so that the system is completely closed and protected until directly prior to use on the patient
- 22F/15M connector on device side as a two-in-one solution
- Standardized connectors for safe connection to other components
- Luer Lock port for gas sampling
- Disposable convenience enhances workflow
- PVC-free design of gas leading components eliminates risk of softeners such as DEHP

A	Catheter mount ErgoStar® CM 40, disposable, 10 cm (3.94 in), double swivel elbow, 50 pcs.	MP01840					
B	Catheter mount ErgoStar® CM 45, disposable, 15 cm (5.91 in), double swivel elbow, 50 pcs.	MP01845					
C	Catheter mount ErgoStar® CM 50, disposable, 10 – 16 cm (3.94 – 6.3 in), extendable, 50 pcs.	MP01850					
D	Catheter mount ErgoStar® CM 55, disposable, 15 cm (5.91 in), straight, 50 pcs.	MP01855					
E	Catheter mount ErgoStar® CM 60, disposable, 10 cm (3.94 in), small, double swivel elbow, 50 pcs.	MP01860					
F	Airway connector ErgoStar® AC 90, disposable, double swivel elbow w/ double cap, 100 pcs.	MP01890					
G	Airway connector ErgoStar® AC 95, disposable, w/ Luer Lock sampling port, 100 pcs.	MP01895					

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



D-12406-2014

B



D-12407-2014

C



D-12406-2014

D



MT-1925-2008

E



D-12409-2014

F



D-12410-2014

G



MT-1926-2008

COMBINATIONS OF FILTER AND TUBE EXTENTION

COMBINATIONS OF FILTER AND TUBE EXTENTION

A	CombiStar Filter HME flex, combination of flexible tube extention MP01850 and filter TwinStar® 55 (MP01805), 25 pcs.	MP04240
B	CombiStar F-HME HEPA flex, combination of flexible tube extention MP01850 and filter TwinStar® HEPA (MP01801), 25 pcs.	MP04242
C	CombiStar mechanical Filter flex, combination of flexible tube extention MP01850 and filter SafeStar® 55 (MP01790), 25 pcs.	MP04246
D	CombiStar A Filter HME, combination of airway connector MP01845 and filter TwinStar® 55 (MP01805), 25 pcs.	MP01671

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



B



C



D



NIV MASKS AND ACCESSORIES

5. NON-INVASIVE PATIENT INTERFACES

NIV MASKS AND ACCESSORIES


ClassicStar® NIV oronasal masks cover the nose and mouth of the patient. Therefore, they provide effective therapy even in case of mouth breathing.

The mask design ensures high patient comfort combined with anatomical fit and effective seal due to individually adjustable mask cushion. The 6-point headgear provides a stable and secure fixation. Adjustable forehead support and forehead pad allow for distributing the pressure from the bridge of nose towards the forehead for more patient comfort. ClassicStar® NIV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NIV-option.

Variants with Anit-Asphyxia Valve (AAV) are intended for Multi-Vendor-Accessory business.

All ClassicStar® NIV masks are for single patient use. Sizing gauges for easy mask size identification as well as a pumpball are delivered together with the masks.

A	NIV oronasal mask ClassicStar® SE, disposable, size S	MP01573					
	NIV oronasal mask ClassicStar® SE, disposable, size M	MP01574					
	NIV oronasal mask ClassicStar® SE, disposable, size L	MP01575					



NIV MASKS AND ACCESSORIES

Oronasal masks are often the first-line strategy in the initial management of acute respiratory failure as they cover nose and mouth. Nasal masks are often chosen if NIV therapy has to be prolonged due to improved comfort and reduced complications associated with oronasal masks. Similar to ClassicStar® oronasal masks, the nasal mask cushion can be individually adjusted to the patient's face providing high patient comfort and effective seal. The 4-point-headgear ensures a stable and secure fixation reducing mask displacement. All ClassicStar® NIV nasal masks are for single patient use. Sizing gauges for easy mask size identification as well as a pumpball are delivered together with the masks.

A NIV nasal mask ClassicStar® non-vented, disposable, size S

MP01623

NIV nasal mask ClassicStar® non-vented, disposable, size M

MP01624

NIV nasal mask ClassicStar® non-vented, disposable, size L

MP01625

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A

D-6409-2009



NIV MASKS AND ACCESSORIES

The ClassicStar® NIV total-face mask seals around the perimeter of the face, covering mouth, nose and eyes. These areas are less pressure sensitive and have smoother facial contours compared to the areas a oronasal masks covers (nasal bridge). The mask is suitable for patients with nasal bridge breakdown or facial irregularities. The air filled cushion can be adjusted individually to the patients face, providing high patient comfort and effective seal.

The 5-point headgear which includes a headstrap at the top of the mask provides a stable and secure attachment, thus preventing mask displacement and reducing pressure. Quick release clips are integrated for quick removal of the mask.

A NIV total-face mask ClassicStar® SE, disposable, size S

MP02910

NIV total-face mask ClassicStar® SE, disposable, size M

MP02911

NIV total-face mask ClassicStar® SE, disposable, size L

MP02912

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



D-6282-2018

NIV MASKS AND ACCESSORIES

The ClassicStar® plus mask series features a soft and anatomically shaped silicone lip as sealing interface to the patients face. The mask has an easy to adjust forehead pad with silicone cushion which allows optimal adaption to the sensitive nasal bridge. The mask is available in four sizes (S, M, L, XL) as oronasal and nasal mask. ClassicStar® plus NIV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NIV-option. Variants with Anti-Asphyxia-Valve (AAV) are intended for Multi-Vendor-Accessory business.

A NIV oronasal mask ClassicStar® plus SE, disposable, size S, 5 pcs.

MP04711

NIV oronasal mask ClassicStar® plus SE, disposable, size M, 5 pcs.

MP04712

NIV oronasal mask ClassicStar® plus SE, disposable, size L, 5 pcs.

MP04713

NIV oronasal mask ClassicStar® plus SE, disposable, size XL, 5 pcs.

MP04714

B NIV nasal mask ClassicStar® plus non vented, disposable, size S, 5 pcs.

MP04741

NIV nasal mask ClassicStar® plus non vented, disposable, size M, 5 pcs.

MP04742

NIV nasal mask ClassicStar® plus non vented, disposable, size L, 5 pcs.

MP04743

NIV nasal mask ClassicStar® plus non vented, disposable, size XL, 5 pcs.

MP04744

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



D-10386-2016

B



D-10389-2016

NIV MASKS AND ACCESSORIES

The NovaStar® plus NiV mask series cover nose and mouth of the patient. Therefore, they provide effective therapy even in case of mouth breathing. This range of products provides ultimate comfort and an effective seal due to the individual adaption of the silicone-gel cushion with a pliable embedded ring. NovaStar® plus NiV masks with standard elbow (SE) have been successfully tested for compatibility with Dräger ventilators which have NiV-option. Variants with Anti-Asphyxia-Valve (AAV) are intended for Multi-Vendor-Accessory business. All NovaStar® plus NiV masks are for multi patient use. Replacement headgears with magnetic clips are available individually (cf. accessories for NiV masks).

A NIV oronasal mask NovaStar® plus, SE, reusable, size S, 3 pcs.

MP05051

NIV oronasal mask NovaStar® plus, SE, reusable, size M, 3 pcs.

MP05052

NIV oronasal mask NovaStar® plus, SE, reusable, size L, 3 pcs.

MP05053

NIV oronasal mask NovaStar® plus, SE, reusable, size XL, 3 pcs.

MP05054

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



D-8729-2019


NIV MASKS AND ACCESSORIES

ACCESSORIES FOR NIV MASKS

A	Two-sided pump ball for air management	MP01590
B	NovaStar® plus Headgear, inclusive magnetic clips, 3 pcs.	MP05080

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A



MT-11178-2006

B



D-25944-2020

ANESTHESIA MASKS AND ACCESSORIES

ANESTHESIA MASKS AND ACCESSORIES

A	① Anesthesia mask LiteStar®, disposable, w/o hook ring, size 1, prenat., 30 pcs.	MP01501
	② Anesthesia mask LiteStar®, disposable, w/o hook ring, size 2, baby, 30 pcs.	MP01502
	③ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 3, infant, 30 pcs.	MP01503
	④ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 4, adult S, 30 pcs.	MP01504
	⑤ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 5, adult M, 30 pcs.	MP01505
	⑥ Anesthesia mask LiteStar®, disposable, w/o hook ring, size 6, adult L, 30 pcs.	MP01506
B	① Hook ring for LiteStar® anesthesia masks, disposable, 15 mm (0.59 in) OD, for sizes 1 – 2, 30 pcs.	MP01551
	② Hook ring for LiteStar® anesthesia masks, disposable, 22 mm (0.87 in) ID, for sizes 3 – 6, 30 pcs.	MP01552

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ANESTHESIA MASKS AND ACCESSORIES

A	① Anesthesia mask ComfortStar®, disposable, hook ring, std., size 0, prenat., 20 pcs.	MP01510
	② Anesthesia mask ComfortStar®, disposable, hook ring, std., size 1, neo., 20 pcs.	MP01511
	③ Anesthesia mask ComfortStar®, disposable, hook ring, std., size 2, baby, 20 pcs.	MP01512
	④ Anesthesia mask ComfortStar®, disposable, hook ring, std., size 3, inf., 20 pcs.	MP01513
	⑤ Anesthesia mask ComfortStar®, disposable, hook ring, std., size 4, adt. S, 20 pcs.	MP01514
	⑥ Anesthesia mask ComfortStar®, disposable, hook ring, std., size 5, adt. M, 20 pcs.	MP01515
	⑦ Anesthesia mask ComfortStar®, disposable, hook ring, std., size 6, adt. L, 20 pcs.	MP01516

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



ANESTHESIA MASKS AND ACCESSORIES

ACCESSORIES FOR ANESTHESIA MASKS

A	① Headgear, disposable, soft fabric, small	MP01559
	② Headgear, disposable, soft fabric, large	MP01560
B	① Headgear, disposable, soft rubber, small	MP01561
	② Headgear, disposable, soft rubber, large	MP01562

EVITA® V800	EVITA® V600	EVITA® INFINITY® V500	EVITA® V300	EVITA® XL	
					A
					
					D-35633-2009
					B
					
					D-35632-2009

HIGH-FLOW THERAPY NASAL CANNULA AND ACCESSORIES

6. HIGH-FLOW THERAPY NASAL CANNULA AND ACCESSORIES


High Flow Therapy allows to provide warm humidified air with enriched oxygen concentration at flow rates up to 60 L per minute to patients via nasal cannulas. Patient comfort is one of the advantages of High Flow therapy over other treatment forms for patient with compromised respiration and can in many cases successfully avoid using less comfortable or even invasive forms of respiratory support.

The HI-Flow Star Nasal Cannula and the inspiratory breathing circuit offer advanced features designed for enhanced patient safety and comfort such as integrated clips in order to provide effective breathing circuit tension relief, a flexible, breathable hose minimizing internal condensation in the high flow nasal cannula as well as soft materials for the prongs and the headgear. Thanks to its symmetrical design, the nasal cannula can be connected from either the left or the right side of the bed. This helps caregivers save time and effort.

The HI-Flow Star products are disposable.

A	① HI-Flow Star Nasal Cannula, disposable, adult S, Flow Range < 60 L/min, 25 pcs.	MP05511	■	■
A	② HI-Flow Star Nasal Cannula, disposable, adult M, Flow Range < 60 L/min, 25 pcs.	MP05512	■	■
A	③ HI-Flow Star Nasal Cannula, disposable, adult L, Flow Range < 60 L/min, 25 pcs.	MP05513	■	■

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL



HIGH-FLOW THERAPY NASAL CANNULA AND ACCESSORIES

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A HI-Flow Star System, heated inspiratory breathing circuit, consists of: heated inspiratory hose, humidifier chamber, breathing hose when used with a ventilator, O ₂ hose when used with an O ₂ blender, bed clip, disposable, 25 pcs.	MP05601
B HI-Flow Star Valve Kit w/ Pop-Up Valve and Connector 22 mm, disposable, 15 pcs.	MP05507
C T-piece adapter for MX300/MX300i and HI-Flow Star System, disposable, 20 pcs.	MP05508



PACK2GO®

7. PACK2GO®

Dräger has developed a comprehensive range of disposable accessories. These accessories can now also be ordered together in a Pack2Go® package, in order to benefit the hospital's costs and logistical/clinical processes. Pack2Go® package solutions include dedicated Dräger and general accessories, which can be used as a system solution with Dräger devices.

A	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, basic + cuvette	MP01656
B	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, basic, 1.8 m	MP01657
C	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, coaxial + cuvette	MP01676
D	Pack2Go® Evita® invasive vent. pack w/ breathing circuit VentStar®, coaxial, 1.8 m	MP01677

EVITA® V800
EVITA® V600
EVITA® INFINITY® V500
EVITA® V300
EVITA® XL

A	
D-98706-2013	
B	
D-98705-2013	
C	
D-10672-2011	
D	
D-10327-2011	

SET2GO

8. SET2GO

With Dräger's Set2go solutions we offer you totally customised sets consisting of an individual selection of single-patient-use accessories. Dräger Set2go makes using and ordering consumables and accessories much simpler for you. You just need one article number to get all the items in a single package and accomplish your specific task.

Please contact your local Dräger sales organisation to create your individual accessory solution.

Not all products, features, or services are for sale in all countries.
Mentioned Trademarks are only registered in certain countries and not necessarily in the country in which this material is released. Go to www.draeger.com/trademarks to find the current status.

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

REGION EUROPE

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION MIDDLE EAST, AFRICA

Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Singapore Pte. Ltd.
61 Science Park Road
The Galen #04-01
Singapore 117525
Tel +65 6872 9288
Fax +65 6259 0398
asia.pacific@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Indústria e Comércio Ltda.
Al. Pucurui - 51 - Tamboré
06406-100 - Barueri - São Paulo
Tel. +55 (11) 4689-4900
relacionamento@draeger.com

CANADA

Draeger Medical Canada Inc.
2425 Skymark Avenue, Unit 1
Mississauga, Ontario, L4W 4Y6
Tel +1 905 212 6600
Toll-free +1 866 343 2273
Fax +1 905 212 6601
Canada.support@draeger.com

Locate your Regional Sales
Representative at:
www.draeger.com/contact

